



Product Specification
Screw cap tube, 15 ml, (LxØ): 120 x 17 mm, PS,
with print

Page 1



Product description

Order number	62.553.542
Product description	Screw cap tube, working volume: 15 ml, (LxØ): 120 x 17 mm, material: PS, conical base, transparent, screw cap, yellow, cap assembled, with print, label/print: white/blue, with graduations, sterile, 50 piece(s)/bag

Product characteristics

Label/ Print	with print
Colour of print/label	white/blue
Graduation	yes
Closure type	screw cap
Base shape	conical base
Volume of work	15 ml
Centrifugation max (RCF)	4100 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Screw cap tube, 15 ml, (LxØ): 120 x 17 mm, PS,
with print**

Page 2

Size

Diameter	17 mm
Length of product	120 mm

Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	yellow

Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	sterile
Sterilisation	Electron irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	24000
Depth of case	425 mm
Width of case	392 mm
Height of case	220 mm
Case volume	0.0367 cbm
Weight of product	0.0085 kg
Weight of case	4.26 kg
EAN of inner box	4038917113007
EAN case	4038917013871

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

