



**Product Specification**  
**Urine Monovette® , 3.2 ml, cap yellow, (LxØ): 75 x 13 mm, 64 piece(s)/bag**

**Page 1**



**Product description**

Order number	10.250.001
Product description	Urine Monovette® , nominal volume: 3.2 ml, without preparation, Urine tube, (LxØ): 75 x 13 mm, transparent, fill mark, material: PP, Luer screw cap, cap yellow, incl. collection tip, with paper label, 64 piece(s)/bag

**Product characteristics**

Type of collection	Aspiration
Type of preparation	without preparation
Label/ Print	with paper label
Colour of print/label	yellow
Graduation	fill mark
Cap	Luer screw cap
Closure type	screw cap
Light protection	no
Application	Urine tube
Base shape	round base
Volume of work	3,2 ml
Recommended centrifugation	400 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Urine Monovette® , 3.2 ml, cap yellow, (LxØ): 75 x 13 mm, 64 piece(s)/bag**

**Page 2**

Centrifugation time 5 Min  
Connection Luer (m)

**Size**

Sample volume 3.2 ml  
Diameter 13 mm  
Length including cap 90 mm  
Length excluding cap 75 mm  
Length including cap and plunger 109 mm

**Material & colours**

Product material Polypropylene (PP)  
Colour of product transparent  
Closure material High Density Polyethylene (HD-PE)  
Colour of cap yellow  
Tube material Polypropylene (PP)  
Piston material High Density Polyethylene (HD-PE)  
Piston plunger material Polystyrene (PS)

**Purity & certification**

Satisfies the requirement IATA, ADR  
Product category In vitrodiagnostic device, CE  
CE certified CE  
Safety data sheet no  
Batched yes

**Packaging**

Minimum order qty. 512  
Type of smallest subpackaging bag  
Piece(s) / inner box 64  
Piece (s) / outer case 512  
Piece(s) / pallet 36864  
Depth of case 388 mm  
Width of case 358 mm  
Height of case 150 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Urine Monovette® , 3.2 ml, cap yellow, (LxØ): 75 x 13 mm, 64 piece(s)/bag**

Case volume	0.0208 cbm
Weight of product	0.0054 kg
Weight of case	3.23 kg
EAN of inner box	4038917279499
EAN case	4038917279482

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

