



Product Specification
Urine Monovette® , 10 ml, cap yellow, (LxØ): 102 x 15 mm, 64 piece(s)/bag

Page 1



Product description

Order number	10.252.001
Product description	Urine Monovette® , nominal volume: 10 ml, without preparation, Urine tube, (LxØ): 102 x 15 mm, transparent, graduated, material: PP, Luer screw cap, cap yellow, incl. collection tip, with paper label, 64 piece(s)/bag

Product characteristics

Type of collection	Aspiration
Type of preparation	without preparation
Label/ Print	with paper label
Colour of print/label	yellow
Graduation	yes
Cap	Luer screw cap
Closure type	screw cap
Light protection	no
Application	Urine tube
Example of use	Urine collection, transport, and storage
Base shape	round base
Volume of work	10 ml

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Urine Monovette® , 10 ml, cap yellow, (LxØ): 102 x 15 mm, 64 piece(s)/bag

Page 2

Recommended centrifugation 400 x g
Centrifugation time 5 Min
Connection Luer (m)

Size

Sample volume 10 ml
Diameter 15 mm
Length including cap 117 mm
Length excluding cap 102 mm
Length including cap and plunger 127 mm

Material & colours

Product material Polypropylene (PP)
Colour of product transparent
Closure material High Density Polyethylene (HD-PE)
Colour of cap yellow
Tube material Polypropylene (PP)
Piston material High Density Polyethylene (HD-PE)
Piston plunger material Polystyrene (PS)

Purity & certification

Satisfies the requirement IATA, ADR
Product category In vitrodiagnostic device, CE
CE certified CE
Batched yes

Packaging

Minimum order qty. 512
Type of smallest bag
subpackaging
Piece(s) / inner box 64
Piece (s) / outer case 512
Piece(s) / pallet 24576
Depth of case 388 mm
Width of case 378 mm
Height of case 215 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Urine Monovette® , 10 ml, cap yellow, (LxØ): 102 x 15 mm, 64 piece(s)/bag

Case volume	0.0315 cbm
Weight of product	0.0085 kg
Weight of case	4.84 kg
EAN of inner box	4038917279932
EAN case	4038917279802

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

