



**Product Specification**  
**Urine Monovette® , Z, cap yellow, (LxØ): 92 x 15 mm, 1 piece(s)/blister**

**Page 1**



**Product description**

Order number	10.258.081
Product description	Urine Monovette® , Z, without preparation, Urine tube, (LxØ): 92 x 15 mm, transparent, graduated, Luer screw cap, cap yellow, incl. collection tip, with paper label, suitable for Tempus® sample transport system, sterile, 1 piece(s)/blister

**Product characteristics**

Type of collection	Aspiration
Type of preparation	without preparation
Label/ Print	with paper label
Colour of print/label	yellow/white
Graduation	yes
Cap	Luer screw cap
Closure type	screw cap
Light protection	no
Application	Urine tube
Suitable for	Tempus® sample transport system
Base shape	round base
Volume of work	8,5 ml

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Urine Monovette® , Z, cap yellow, (LxØ): 92 x 15 mm, 1 piece(s)/blister**

**Page 2**

Recommended centrifugation 400 x g  
Centrifugation time 5 Min  
Connection Luer (m)

**Size**

Diameter 15 mm  
Length including cap 106 mm  
Length excluding cap 92 mm  
Length including cap and plunger 127 mm

**Material & colours**

Colour of product transparent  
Closure material High Density Polyethylene (HD-PE)  
Colour of cap yellow  
Tube material Polypropylene (PP)  
Piston material High Density Polyethylene (HD-PE)  
Piston plunger material Polystyrene (PS)

**Purity & certification**

Satisfies the requirement IATA, ADR  
Product category In-vitro diagnostic (IVD)  
CE certified CE - manufacturer's self-declaration  
Purity standard sterile  
Sterilisation Electron irradiation  
Safety data sheet no  
Batched yes

**Packaging**

Minimum order qty. 500  
Type of smallest blister  
subpackaging  
Piece(s) / inner box 100  
Piece (s) / outer case 500  
Piece(s) / pallet 13500  
Depth of box 351 mm  
Width of box 151 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Urine Monovette® , Z, cap yellow, (LxØ): 92 x 15 mm, 1 piece(s)/blister**

**Page 3**

Height of box	171 mm
Depth of case	778 mm
Width of case	368 mm
Height of case	195 mm
Case volume	0.0558 cbm
Weight of product	0.0078 kg
Weight of case	5.51 kg
EAN of inner box	4038917992558
EAN case	4038917992565

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2025-12-13

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**