



**Product description**

Order number	11.1240.100
Product description	Urine transfer unit for V-Monovette® Urine, for Urine transfer, 50 piece(s)/minigrip bag

**Product characteristics**

Type of collection	Vacuum
Light protection	no
Design	With wings on the tube holder
Application	Urine transfer
Suitable for	V-Monovette® Urine
Connection	Vacuum tube

**Size**

Length of tube	10 cm
----------------	-------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Colour of product yellow

### **Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE supplier
Safety data sheet	no
Batched	yes

### **Packaging**

Minimum order qty.	500
Type of smallest subpackaging	minigrip bag
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	30000
Depth of box	300 mm
Width of box	250 mm
Height of box	10 mm
Depth of case	380 mm
Width of case	350 mm
Height of case	170 mm
Case volume	0.0226 cbm
Weight of product	0.002 kg
Weight of case	1.8 kg
EAN of inner box	6947462401157
EAN case	6947462401423

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

