



## Product description

Order number	11.2253.001
Product description	V-Monovette® Urine, Boric acid, nominal volume: 4 ml, Boric acid, vacuum system for closed urine transfer, (LxØ): 75 x 13 mm, transparent, fill mark, material: PET, membrane cap, cap green, with paper label, sterile, 50 piece(s)/bag

## Product characteristics

Type of collection	Vacuum
Type of preparation	Boric acid
Preparation concentration	15 mg/ml
Label/ Print	with paper label
Colour of print/label	green
Graduation	fill mark
Cap	membrane cap
Closure type	push cap
Light protection	no
Application	vacuum system for closed urine transfer
Base shape	round base

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Product Specification

**V-Monovette® Urine, Boric acid, 4 ml, cap green,  
(LxØ): 75 x 13 mm, 50 piece(s)/bag**

**Page 2**

Volume of work 4 ml  
Connection V-Monovette®

### Size

Sample volume 4 ml  
Diameter 13 mm  
Length including cap 81 mm  
Length excluding cap 75 mm

### Material & colours

Product material Polyethylene therephthalate (PET)  
Colour of product transparent  
Closure material High Density Polyethylene (HD-PE)  
Colour of cap green  
Membrane material rubber  
Tube material Polyethylene therephthalate (PET)

### Purity & certification

Product category In-vitro diagnostic (IVD)  
CE certified CE - manufacturer's self-declaration  
Purity standard sterile  
Sterilisation Gamma irradiation  
Safety data sheet yes  
Batched yes

### Packaging

Minimum order qty. 500  
Type of smallest bag  
subpackaging  
Piece(s) / inner box 50  
Piece (s) / outer case 500  
Piece(s) / pallet 39000  
Depth of case 398 mm  
Width of case 368 mm  
Height of case 130 mm  
Case volume 0.019 cbm  
Weight of product 0.006 kg

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2025-12-13

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**

**V-Monovette® Urine, Boric acid, 4 ml, cap green,  
(LxØ): 75 x 13 mm, 50 piece(s)/bag**

Weight of case	3.47 kg
EAN of inner box	4038917245128
EAN case	4038917245111

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

