



**Product Specification**  
**V-Monovette® Urine, 9.5 ml, cap yellow, (LxØ):**  
**100 x 15 mm, 50 piece(s)/bag**

**Page 1**



**Product description**

Order number	11.2552.001
Product description	V-Monovette® Urine, nominal volume: 9.5 ml, without preparation, vacuum system for closed urine transfer, (LxØ): 100 x 15 mm, transparent, fill mark, material: PET, membrane cap, cap yellow, with paper label, sterile, 50 piece(s)/bag

**Product characteristics**

Type of collection	Vacuum
Type of preparation	without preparation
Label/ Print	with paper label
Colour of print/label	yellow
Graduation	fill mark
Cap	membrane cap
Closure type	push cap
Light protection	no
Application	vacuum system for closed urine transfer
Base shape	conical base
Volume of work	9,5 ml

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**V-Monovette® Urine, 9.5 ml, cap yellow, (LxØ):**  
**100 x 15 mm, 50 piece(s)/bag**

**Page 2**

Recommended centrifugation 400 x g  
Centrifugation time 5 Min

**Size**

Sample volume 9.5 ml  
Diameter 15 mm  
Length including cap 106 mm  
Length excluding cap 100 mm

**Material & colours**

Product material Polyethylene therephthalate (PET)  
Colour of product transparent  
Closure material High Density Polyethylene (HD-PE)  
Colour of cap yellow  
Membrane material rubber  
Tube material Polyethylene therephthalate (PET)

**Purity & certification**

Product category In-vitro diagnostic (IVD)  
CE certified CE - manufacturer's self-declaration  
Purity standard sterile  
Sterilisation Gamma irradiation  
Safety data sheet no  
Batched yes

**Packaging**

Minimum order qty. 500  
Type of smallest bag  
subpackaging  
Piece(s) / inner box 50  
Piece (s) / outer case 500  
Piece(s) / pallet 22000  
Depth of case 503 mm  
Width of case 375 mm  
Height of case 155 mm  
Case volume 0.0292 cbm  
Weight of product 0.0088 kg

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**V-Monovette® Urine, 9.5 ml, cap yellow, (LxØ):**  
**100 x 15 mm, 50 piece(s)/bag**

**Page 3**

Weight of case	4.88 kg
EAN of inner box	4038917244886
EAN case	4038917244879

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2025-12-14

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**