




Product Specification
Seraplas® valve filter, white, for separation of serum/plasma from the blood cells after centrifugation

Page 1



Product description

Order number	53.677
Product description	Seraplas® valve filter, V11, white, for separation of serum/plasma from the blood cells after centrifugation, 500 piece(s)/bag

Product characteristics

Application	venous blood collection
Suitable for	separation of serum/plasma from the blood cells after centrifugation

Size

Length of product	92 mm
-------------------	-------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Seraplas® valve filter, white, for separation of
serum/plasma from the blood cells after
centrifugation**

Page 2

Material & colours

Product material	Polyethylene (PE)
Colour of product	white

Purity & certification

Product category	no medical device IVD
Safety data sheet	no
Batched	no

Packaging

Minimum order qty.	5000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	5000
Piece(s) / pallet	120000
Depth of case	398 mm
Width of case	398 mm
Height of case	415 mm
Case volume	0.0657 cbm
Weight of product	0.0015 kg
Weight of case	8.5 kg
EAN of inner box	4038917083317
EAN case	4038917002257

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-27

This document was prepared by EDP
support and is valid without signature.



SARSTEDT