



Product description

Order number	55.1570
Product description	Adapter tube, (LxØ): 54 x 11 mm, material: PP, transparent, 250 piece(s)/bag

Product characteristics

Suitable for	S-Monovette® with Ø 8 mm
Base shape	flat base

Size

Diameter	11 mm
Length of product	54 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 1000
Type of smallest
subpackaging bag
Piece(s) / inner box 250
Piece (s) / outer case 1000
Piece(s) / pallet 78000
Depth of case 428 mm
Width of case 323 mm
Height of case 130 mm
Case volume 0.018 cbm
Weight of product 0.0011 kg
Weight of case 1.51 kg
EAN of inner box 4038917087063
EAN case 4038917005531

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

