



Product description

| | |
|---------------------|---|
| Order number | 55.1571.051 |
| Product description | Adapter tube, (LxØ): 65 x 13 mm, material: PP, light blue, 250 piece(s)/bag |

Product characteristics

| | |
|--------------|---|
| Suitable for | S-Monovette® with Ø 8 mm, Microvette® and micro tubes |
| Base shape | round base |

Size

| | |
|-------------------|-------|
| Diameter | 13 mm |
| Length of product | 65 mm |

Material & colours

| | |
|-------------------|--------------------|
| Product material | Polypropylene (PP) |
| Colour of product | light blue |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| Batched | no |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1000 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 250 |
| Piece (s) / outer case | 1000 |
| Piece(s) / pallet | 72000 |
| Depth of case | 333 mm |
| Width of case | 303 mm |
| Height of case | 190 mm |
| Case volume | 0.0192 cbm |
| Weight of product | 0.0019 kg |
| Weight of case | 2.23 kg |
| EAN of inner box | 4038917283380 |
| EAN case | 4038917283373 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

