



### **Product description**

Order number	55.1571
Product description	Adapter tube, (LxØ): 65 x 13 mm, material: PP, transparent, 250 piece(s)/bag

### **Product characteristics**

Closure type	screw cap
Suitable for	S-Monovette® with Ø 8 mm, Microvette® and micro tubes
Base shape	round base

### **Size**

Diameter	13 mm
Length of product	65 mm

### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Purity & certification**

Product category	no medical device   IVD
Batched	no

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	250
Piece (s) / outer case	1000
Piece(s) / pallet	72000
Depth of case	333 mm
Width of case	303 mm
Height of case	190 mm
Case volume	0.0192 cbm
Weight of product	0.0019 kg
Weight of case	2.23 kg
EAN of inner box	4038917242523
EAN case	4038917242516

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

