



### **Product description**

Order number	55.1577.002
Product description	Carrier tube, (LxØ): 60 x 11.5 mm, material: PP, brown, 250 piece(s)/bag

### **Product characteristics**

Light protection	yes
Suitable for	Multivette®
Base shape	flat base

### **Size**

Diameter	11.5 mm
Length of product	60 mm

### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	brown

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Purity & certification**

Product category	no medical device   IVD
Batched	no

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	250
Piece (s) / outer case	1000
Piece(s) / pallet	78000
Depth of case	428 mm
Width of case	323 mm
Height of case	130 mm
Case volume	0.018 cbm
Weight of product	0.0013 kg
Weight of case	1.7 kg
EAN of inner box	4038917131032
EAN case	4038917023597

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

