



Product description

Order number	55.1577
Product description	Carrier tube, (LxØ): 60 x 11.5 mm, material: PP, transparent, 250 piece(s)/bag

Product characteristics

Closure type	screw cap
Suitable for	Multivette®
Base shape	flat base

Size

Diameter	11.5 mm
Length of product	60 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 1000
Type of smallest bag
subpackaging
Piece(s) / inner box 250
Piece (s) / outer case 1000
Piece(s) / pallet 78000
Depth of case 428 mm
Width of case 323 mm
Height of case 130 mm
Case volume 0.018 cbm
Weight of product 0.0013 kg
Weight of case 1.7 kg
EAN of inner box 4038917087070
EAN case 4038917005548

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

