



### **Product description**

Order number	55.525.003
Product description	Tube, working volume: 5 ml, (LxØ): 75 x 13 mm, material: PP, round base, transparent, suitable for Microvette® and micro tube with screw cap, 250 piece(s)/bag

### **Product characteristics**

Suitable for	Microvette® and micro tube with screw cap
Base shape	round base
Volume of work	5 ml

### **Size**

Diameter	13 mm
Length of product	75 mm

### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Tube, 5 ml, (LxØ): 75 x 13 mm, PP**

**Page 2**

**Purity & certification**

Product category	no medical device   IVD
Batched	no

**Packaging**

Minimum order qty.	2000
Type of smallest subpackaging	bag
Piece(s) / inner box	250
Piece (s) / outer case	2000
Piece(s) / pallet	72000
Depth of case	398 mm
Width of case	398 mm
Height of case	265 mm
Case volume	0.042 cbm
Weight of product	0.0026 kg
Weight of case	5.7 kg
EAN of inner box	4038917109154
EAN case	4038917012362

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-29

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**