



Product description

Order number	55.526.053
Product description	Tube, working volume: 5 ml, (LxØ): 75 x 12 mm, material: PP, round base, transparent, push cap, suitable for Abbott m1000, 500 piece(s)/StackPack

Product characteristics

Graduation	no
Closure type	push cap
Suitable for	Abbott m1000
Base shape	round base
Volume of work	5 ml
Centrifugation max (RCF)	5000 x g

Size

Diameter	12 mm
Length of product	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Tube, 5 ml, (LxØ): 75 x 12 mm, PP

Page 2

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	2000
Type of smallest subpackaging	StackPack
Piece(s) / inner box	500
Piece (s) / outer case	2000
Piece(s) / pallet	48000
Depth of case	378 mm
Width of case	210 mm
Height of case	391 mm
Case volume	0.0309 cbm
Weight of product	0.0026 kg
Weight of case	6 kg
EAN of inner box	4038917109710
EAN case	4038917080033

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-15

This document was prepared by EDP support and is valid without signature.



SARSTEDT