



### **Product description**

Order number 65.1379.002  
Product description Cap, red, suitable for tubes Ø 12-17 mm, 1,000 piece(s)/bag

### **Product characteristics**

Cap standard cap  
Suitable for tubes Ø 12-17 mm

### **Material & colours**

Product material Low Density Polyethylene (LD-PE)  
Colour of product red

### **Purity & certification**

Product category no medical device | IVD  
Batched yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Packaging**

Minimum order qty.	10000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	10000
Piece(s) / pallet	240000
Depth of case	413 mm
Width of case	257 mm
Height of case	464 mm
Case volume	0.0492 cbm
Weight of product	0.0004 kg
Weight of case	5.36 kg
EAN of inner box	4038917127981
EAN case	4038917021517

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

