



**Product Specification**  
**Insert, suitable for tubes and S-Monovette® Ø 15 mm, transparent**

**Page 1**



**Product description**

Order number	73.1052
Product description	Insert, suitable for tubes and S-Monovette® Ø 15 mm, material: PS, length: 23.5 mm, transparent, 500 piece(s)/bag

**Product characteristics**

Graduation	no
Suitable for	tubes and S-Monovette® Ø 15 mm
Base shape	skirted conical base

**Size**

Diameter of opening	11 mm
Length of product	23.5 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Insert, suitable for tubes and S-Monovette® Ø 15**  
**mm, transparent**

**Page 2**

**Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent

**Purity & certification**

Product category	no medical device   IVD
Batched	no

**Packaging**

Minimum order qty.	2000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	2000
Piece(s) / pallet	156000
Depth of case	428 mm
Width of case	323 mm
Height of case	130 mm
Case volume	0.018 cbm
Weight of product	0.0008 kg
Weight of case	2.02 kg
EAN of inner box	4038917085151
EAN case	4038917003834

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-03

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**