




Product Specification
Wall holder, for tournistrip® disposable
tourniquet (item no 95.1006)

Page 1



Product description

Order number	95.1005
Product description	Wall holder, for tournistrip® disposable tourniquet (item no 95.1006)

Product characteristics

Suitable for	tournistrip® disposable tourniquet (item no 95.1006)
--------------	--

Material & colours

Product material	metal
------------------	-------

Purity & certification

Product category	no medical device IVD
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Wall holder, for tournistrip® disposable
tourniquet (item no 95.1006)

Page 2

Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	324
Depth of case	268 mm
Width of case	191 mm
Height of case	100 mm
Case volume	0.0051 cbm
Weight of product	0.21 kg
Weight of case	0.33 kg
EAN of inner box	4038917348560
EAN case	4038917348553

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-20

This document was prepared by EDP support and is valid without signature.



SARSTEDT