



Product description

Order number	95.1007.101
Product description	Tourniquet, latex-free

Product characteristics

Application	venous blood collection
-------------	-------------------------

Purity & certification

Product category	no medical device   IVD
Latex-free	yes
Batched	yes

Packaging

Minimum order qty.	250
Type of smallest subpackaging	case
Piece(s) / inner box	250
Piece (s) / outer case	2000
Piece(s) / pallet	42000
Case volume	0 cbm
Weight of case	7.01 kg

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Tourniquet, latex-free**

**Page 2**

EAN of inner box

4038917263214

EAN case

4038917398220

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**