



### **Product description**

Order number	14.1205.100
Product description	Multi-adapter, Luer lock, connection: Luer lock, for S-Monovette®, 1 piece(s)/blister, 100 piece(s)/case, sterile, non-pyrogenic/endotoxin-free

### **Product characteristics**

Design	Luer lock
Suitable for	S-Monovette®
Connection	Luer lock

### **Size**

Length of product	47 mm
-------------------	-------

### **Material & colours**

Product material	Polyethylene (PE)/acrylonitrile butadiene styrene (ABS)
Colour of product	translucent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Purity & certification**

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	sterile, non-pyrogenic/endotoxin-free
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

### **Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	blister
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	80000
Depth of box	268 mm
Width of box	81 mm
Height of box	71 mm
Depth of case	373 mm
Width of case	287 mm
Height of case	172 mm
Case volume	0.0184 cbm
Weight of product	0.0021 kg
Weight of case	3.11 kg
EAN of inner box	4038917126700
EAN case	4038917020817

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

