



### **Product description**

|                     |   |
|---------------------|---|
| Order number        | 14.1205.100   |
| Product description | Multi-adapter, Luer lock, connection: Luer lock, for S-Monovette®, 1 piece(s)/blister, 100 piece(s)/case, sterile, non-pyrogenic/endotoxin-free |

### **Product characteristics**

|              |              |
|--------------|--------------|
| Design       | Luer lock    |
| Suitable for | S-Monovette® |
| Connection   | Luer lock    |

### **Size**

|                   |       |
|-------------------|-------|
| Length of product | 47 mm |
|-------------------|-------|

### **Material & colours**

|                   |   |
|-------------------|---|
| Product material  | Polyethylene (PE)/acrylonitrile butadiene styrene (ABS) |
| Colour of product | translucent   |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Purity & certification**

|                  |                                       |
|------------------|---------------------------------------|
| Product category | In vitro diagnostic device, CE        |
| CE certified     | CE                                    |
| Purity standard  | sterile, non-pyrogenic/endotoxin-free |
| Sterilisation    | Gamma irradiation                     |
| Pyrogen-free     | yes                                   |
| Batched          | yes                                   |

### **Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 1000          |
| Type of smallest subpackaging | blister       |
| Piece(s) / inner box          | 100           |
| Piece (s) / outer case        | 1000          |
| Piece(s) / pallet             | 80000         |
| Depth of box                  | 268 mm        |
| Width of box                  | 81 mm         |
| Height of box                 | 71 mm         |
| Depth of case                 | 373 mm        |
| Width of case                 | 178 mm        |
| Height of case                | 287 mm        |
| Case volume                   | 0.0191 cbm    |
| Weight of product             | 0.0021 kg     |
| Weight of case                | 3.11 kg       |
| EAN of inner box              | 4038917126700 |
| EAN case                      | 4038917020817 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

