



Product description

| | |
|---------------------|---|
| Order number | 14.1207 |
| Product description | Blood Culture Adapter LongNeck, For narrow bottle neck, membrane screw cap assembled, connection: blood culture bottle, for use with (Safety) Multifly® needle, 1 piece(s)/blister, 50 piece(s)/case, sterile, non-pyrogenic/endotoxin-free |

Product characteristics

| | |
|--------------|------------------------------------|
| Design | membrane screw cap assembled |
| Suitable for | use with (Safety) Multifly® needle |
| Connection | blood culture bottle |

Size

| | |
|---------------------|-------|
| Diameter of opening | 19 mm |
| Length of product | 74 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Colour of product transparent

Purity & certification

Product category In vitrodiagnostic device, CE
CE certified CE
Purity standard sterile, non-pyrogenic/endotoxin-free
Sterilisation Electron irradiation
Batched yes

Packaging

Minimum order qty. 250
Type of smallest blister
subpackaging
Piece(s) / inner box 50
Piece (s) / outer case 250
Piece(s) / pallet 9000
Depth of case 578 mm
Width of case 368 mm
Height of case 195 mm
Case volume 0.0415 cbm
Weight of product 0.0041 kg
Weight of case 2.14 kg
EAN of inner box 4038917364188
EAN case 4038917364171

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

