



**Product Specification**  
**Universal blood culture adapter, For wide or narrow bottle neck, membrane adapter assembled**

**Page 1**



**Product description**

Order number	14.1209
Product description	Universal blood culture adapter, For wide or narrow bottle neck, membrane adapter assembled, connection: blood culture bottle, for use with (Safety) Multifly® needle, 1 piece(s)/blister, 20 piece(s)/case, sterile, non-pyrogenic/endotoxin-free

**Product characteristics**

Design	membrane adapter assembled
Suitable for	use with (Safety) Multifly® needle
Connection	blood culture bottle

**Size**

Diameter of opening	38 mm
Length of product	84 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Universal blood culture adapter, For wide or narrow bottle neck, membrane adapter assembled**

**Page 2**

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

**Purity & certification**

Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	sterile, non-pyrogenic/endotoxin-free
Sterilisation	Electron irradiation
Batched	yes

**Packaging**

Minimum order qty.	100
Type of smallest subpackaging	blister
Piece(s) / inner box	20
Piece (s) / outer case	100
Piece(s) / pallet	3600
Depth of case	578 mm
Width of case	368 mm
Height of case	195 mm
Case volume	0.0415 cbm
Weight of product	0.0079 kg
Weight of case	1.94 kg
EAN of inner box	4038917354776
EAN case	4038917354769

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-04

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**