



**Product Specification**  
**S-Monovette® needle, 21G x 1 1/2", green, 1  
piece(s)/blister**

**Page 1**



**Product description**

|                     |   |
|---------------------|---|
| Order number        | 85.1162   |
| Product description | S-Monovette® needle, needle (ØxL): 21G x 1 1/2", coding colour: green, 1 piece(s)/blister, sterile, non-pyrogenic/endotoxin-free, 100 piece(s)/case |

**Product characteristics**

|                |                    |
|----------------|--------------------|
| Coding colour  | green              |
| Safety product | no                 |
| Connection     | membrane screw cap |

**Size**

|                               |               |
|-------------------------------|---------------|
| Diameter of needle (external) | 21 G / 0,8 mm |
| Length of needle              | 38 mm         |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**S-Monovette® needle, 21G x 1 1/2", green, 1  
piece(s)/blister**

**Page 2**

**Purity & certification**

|                  |                                       |
|------------------|---------------------------------------|
| Product category | Medical device                        |
| CE certified     | CE 0124 (DEKRA)                       |
| Purity standard  | sterile, non-pyrogenic/endotoxin-free |
| Sterilisation    | Gamma irradiation                     |
| Pyrogen-free     | yes                                   |
| Batched          | yes                                   |

**Packaging**

|                                  |               |
|----------------------------------|---------------|
| Minimum order qty.               | 1000          |
| Type of smallest<br>subpackaging | blister       |
| Piece(s) / inner box             | 100           |
| Piece (s) / outer case           | 1000          |
| Piece(s) / pallet                | 78000         |
| Depth of box                     | 268 mm        |
| Width of box                     | 111 mm        |
| Height of box                    | 71 mm         |
| Depth of case                    | 373 mm        |
| Width of case                    | 233 mm        |
| Height of case                   | 290 mm        |
| Case volume                      | 0.0252 cbm    |
| Weight of product                | 0.0017 kg     |
| Weight of case                   | 2.89 kg       |
| EAN of inner box                 | 4038917085526 |
| EAN case                         | 4038917004176 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-08

This document was prepared by EDP  
support and is valid without signature.



**SARSTEDT**