



## Product description

|                     |                                                                                                                                                             |
|---------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Order number        | 85.1441.200                                                                                                                                                 |
| Product description | S-Monovette® safety needle, needle (ØxL): 22G x 1",<br>coding colour: black, 1 piece(s)/blister, sterile,<br>non-pyrogenic/endotoxin-free, 50 piece(s)/case |

## Product characteristics

|                |                    |
|----------------|--------------------|
| Coding colour  | black              |
| Safety product | yes                |
| Connection     | membrane screw cap |

## Size

|                               |               |
|-------------------------------|---------------|
| Diameter of needle (external) | 22 G / 0,7 mm |
| Length of needle              | 25.4 mm       |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Product Specification

**S-Monovette® safety needle, 22G x 1", black, 1 piece(s)/blister**

**Page 2**

### Purity & certification

|                  |                                       |
|------------------|---------------------------------------|
| Product category | Medical device (MD)                   |
| CE certified     | CE0124 - DEKRA                        |
| Purity standard  | sterile, non-pyrogenic/endotoxin-free |
| Sterilisation    | Gamma irradiation                     |
| Pyrogen-free     | yes                                   |
| Batched          | yes                                   |

### Packaging

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 500           |
| Type of smallest subpackaging | blister       |
| Piece(s) / inner box          | 50            |
| Piece (s) / outer case        | 500           |
| Piece(s) / pallet             | 39000         |
| Depth of box                  | 268 mm        |
| Width of box                  | 111 mm        |
| Height of box                 | 71 mm         |
| Depth of case                 | 373 mm        |
| Width of case                 | 233 mm        |
| Height of case                | 290 mm        |
| Case volume                   | 0.0252 cbm    |
| Weight of product             | 0.0026 kg     |
| Weight of case                | 2.46 kg       |
| EAN of inner box              | 4038917129206 |
| EAN case                      | 4038917361972 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2025-12-14

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**