



## Product Specification

Sample tube, Citrate 9NC 0.106 mol/l 3.2%, 5 ml,  
cap green, (LxØ): 75 x 13 mm, with print

Page 1



### Product description

|                     |   |
|---------------------|---|
| Order number        | 32.525.326  |
| Product description | Sample tube, nominal volume: 5 ml, cap green, preparation: Citrate 9NC, push cap, (LxØ): 75 x 13 mm, round base, with print, 100 piece(s)/polystyrene container |

### Product characteristics

|                            |             |
|----------------------------|-------------|
| Type of preparation        | Citrate 9NC |
| Preparation concentration  | 0.106 mol/l |
| Colour code                | EU          |
| Label/ Print               | with print  |
| Colour of print/label      | green       |
| Graduation                 | fill mark   |
| Closure type               | push cap    |
| Base shape                 | round base  |
| Recommended centrifugation | 1800 x g    |
| Centrifugation time        | 10 Min      |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Product Specification

**Sample tube, Citrate 9NC 0.106 mol/l 3.2%, 5 ml,  
cap green, (LxØ): 75 x 13 mm, with print**

**Page 2**

### Size

|                      |       |
|----------------------|-------|
| Sample volume        | 5 ml  |
| Diameter             | 13 mm |
| Length of product    | 75 mm |
| Length excluding cap | 75 mm |

### Material & colours

|                   |                                  |
|-------------------|----------------------------------|
| Colour of product | transparent                      |
| Closure material  | Low Density Polyethylene (LD-PE) |
| Colour of cap     | green                            |
| Tube material     | Polypropylene (PP)               |

### Purity & certification

|                  |                                |
|------------------|--------------------------------|
| Product category | In vitro diagnostic device, CE |
| CE certified     | CE                             |
| Batched          | yes                            |

### Packaging

|                               |                       |
|-------------------------------|-----------------------|
| Minimum order qty.            | 500                   |
| Type of smallest subpackaging | polystyrene container |
| Piece(s) / inner box          | 100                   |
| Piece (s) / outer case        | 500                   |
| Piece(s) / pallet             | 37500                 |
| Depth of case                 | 458 mm                |
| Width of case                 | 348 mm                |
| Height of case                | 115 mm                |
| Case volume                   | 0.0183 cbm            |
| Weight of product             | 0.0038 kg             |
| Weight of case                | 2.42 kg               |
| EAN of inner box              | 4038917109277         |
| EAN case                      | 4038917012454         |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-03

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**