



**Product Specification**  
**Sample tube, EDTA K3E, 4 ml, cap violet, (LxØ):  
75 x 12 mm, with paper label**

**Page 1**



**Product description**

Order number	36.338.100
Product description	Sample tube, nominal volume: 4 ml, cap violet, preparation: K3 EDTA, push cap, (LxØ): 75 x 12 mm, round base, with paper label, 100 piece(s)/bag

**Product characteristics**

Type of preparation	K3 EDTA
Preparation concentration	1.6 mg/ml blood
Colour code	EU
Label/ Print	with paper label
Colour of print/label	violet
Graduation	fill mark
Cap	Skirted push cap
Closure type	push cap
Base shape	round base

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Sample tube, EDTA K3E, 4 ml, cap violet, (LxØ):**  
**75 x 12 mm, with paper label**

**Page 2**

**Size**

Sample volume	4 ml
Diameter	12 mm
Length of product	75 mm
Length excluding cap	75 mm

**Material & colours**

Colour of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Colour of cap	violet
Tube material	Polypropylene (PP)

**Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

**Packaging**

Minimum order qty.	50000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	2000
Piece(s) / pallet	60000
Depth of case	398 mm
Width of case	398 mm
Height of case	340 mm
Case volume	0.0539 cbm
Weight of product	0.0034 kg
Weight of case	7.48 kg
EAN of inner box	4038917105545
EAN case	4038917010658

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

