




Product Specification

Sample tube, Fluoride/heparin FH, 2 ml, cap yellow, (LxØ): 75 x 12 mm, with print

Page 1



Product description

| | |
|---------------------|---|
| Order number | 36.339 |
| Product description | Sample tube, nominal volume: 2 ml, cap yellow, preparation: Fluoride + heparin, push cap, (LxØ): 75 x 12 mm, round base, with print, 100 piece(s)/bag |

Product characteristics

| | |
|----------------------------|----------------------------------|
| Type of preparation | Fluoride + heparin |
| Preparation concentration | 1.0 mg/ml blood + 16 IU/ml blood |
| Colour code | EU |
| Label/ Print | with print |
| Colour of print/label | yellow |
| Graduation | yes |
| Cap | Skirted push cap |
| Closure type | push cap |
| Base shape | round base |
| Recommended centrifugation | 2000 x g |
| Centrifugation time | 10 Min |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Sample tube, Fluoride/heparin FH, 2 ml, cap yellow, (LxØ): 75 x 12 mm, with print

Page 2

Size

| | |
|----------------------|-------|
| Sample volume | 2 ml |
| Diameter | 12 mm |
| Length of product | 75 mm |
| Length excluding cap | 75 mm |

Material & colours

| | |
|-------------------|----------------------------------|
| Colour of product | transparent |
| Closure material | Low Density Polyethylene (LD-PE) |
| Colour of cap | yellow |
| Tube material | Polypropylene (PP) |

Purity & certification

| | |
|------------------|--------------------------------------|
| Product category | In-vitro diagnostic (IVD) |
| CE certified | CE - manufacturer's self-declaration |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 2000 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 100 |
| Piece (s) / outer case | 2000 |
| Piece(s) / pallet | 60000 |
| Depth of case | 398 mm |
| Width of case | 398 mm |
| Height of case | 340 mm |
| Case volume | 0.0539 cbm |
| Weight of product | 0.0033 kg |
| Weight of case | 7.31 kg |
| EAN of inner box | 4038917081245 |
| EAN case | 4038917000383 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

