




**Product Specification**  
**Sample tube, Fluoride/heparin FH, 2 ml, cap yellow, (LxØ): 75 x 12 mm, with print**

**Page 1**



**Product description**

|                     |                                                                                                                                                       |
|---------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------|
| Order number        | 36.339                                                                                                                                                |
| Product description | Sample tube, nominal volume: 2 ml, cap yellow, preparation: Fluoride + heparin, push cap, (LxØ): 75 x 12 mm, round base, with print, 100 piece(s)/bag |

**Product characteristics**

|                            |                                  |
|----------------------------|----------------------------------|
| Type of preparation        | Fluoride + heparin               |
| Preparation concentration  | 1.0 mg/ml blood + 16 IU/ml blood |
| Colour code                | EU                               |
| Label/ Print               | with print                       |
| Colour of print/label      | yellow                           |
| Graduation                 | yes                              |
| Cap                        | Skirted push cap                 |
| Closure type               | push cap                         |
| Base shape                 | round base                       |
| Recommended centrifugation | 2000 x g                         |
| Centrifugation time        | 10 Min                           |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Sample tube, Fluoride/heparin FH, 2 ml, cap yellow, (LxØ): 75 x 12 mm, with print**

**Page 2**

**Size**

|                      |       |
|----------------------|-------|
| Sample volume        | 2 ml  |
| Diameter             | 12 mm |
| Length of product    | 75 mm |
| Length excluding cap | 75 mm |

**Material & colours**

|                   |                                  |
|-------------------|----------------------------------|
| Colour of product | transparent                      |
| Closure material  | Low Density Polyethylene (LD-PE) |
| Colour of cap     | yellow                           |
| Tube material     | Polypropylene (PP)               |

**Purity & certification**

|                  |                               |
|------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified     | CE                            |
| Batched          | yes                           |

**Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 2000          |
| Type of smallest subpackaging | bag           |
| Piece(s) / inner box          | 100           |
| Piece (s) / outer case        | 2000          |
| Piece(s) / pallet             | 60000         |
| Depth of case                 | 398 mm        |
| Width of case                 | 398 mm        |
| Height of case                | 340 mm        |
| Case volume                   | 0.0539 cbm    |
| Weight of product             | 0.0033 kg     |
| Weight of case                | 7.31 kg       |
| EAN of inner box              | 4038917081245 |
| EAN case                      | 4038917000383 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

