



## Product description

Order number	41.1350.005
Product description	Micro sample tube Citrate 9NC 0.106 mol/l 3.2%, venous blood collection, preparation: Citrate 9NC, nominal volume: 1.3 ml, (LxØ) with cap: 47 x 10.8 mm, screw cap, cap: green, colour code EU, with paper label, label/print: green, skirted conical base, 100 piece(s)/bag

## Product characteristics

Type of collection	venous
Type of preparation	Citrate 9NC
Preparation concentration	0.106 mol/l
Colour code	EU
Label/ Print	with paper label
Colour of print/label	green
Graduation	yes
Closure type	screw cap
Application	venous blood collection
Base shape	skirted conical base
Volume of work	1,3 ml

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Micro sample tube Citrate 9NC 0.106 mol/l 3.2%,  
1.3 ml, screw cap, EU**

**Page 2**

Recommended centrifugation 1500 U/min  
Centrifugation time 10 Min

**Size**

Sample volume 1.3 ml  
Diameter 10.8 mm  
Length of product 47 mm  
Length including cap 47 mm  
Length excluding cap 45 mm

**Material & colours**

Colour of product transparent  
Closure material Polypropylene (PP)  
Colour of cap green  
Tube material Polypropylene (PP)

**Purity & certification**

Satisfies the requirement IATA, ADR  
Product category In vitro diagnostic device, CE  
CE certified CE  
Batched yes

**Packaging**

Minimum order qty. 1000  
Type of smallest bag  
subpackaging  
Piece(s) / inner box 100  
Piece (s) / outer case 1000  
Piece(s) / pallet 96000  
Depth of case 348 mm  
Width of case 278 mm  
Height of case 140 mm  
Case volume 0.0135 cbm  
Weight of product 0.0017 kg  
Weight of case 2.02 kg  
EAN of inner box 4038917127899  
EAN case 4038917021463

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

