



Product description

Order number	41.1500.005
Product description	Micro sample tube Serum Gel CAT, venous blood collection, preparation: Clotting activator/gel, nominal volume: 1.1 ml, (LxØ) with cap: 44 x 10.8 mm, push cap, cap: brown, colour code EU/ISO, with plastic label, label/print: brown, skirted conical base, 100 piece(s)/bag

Product characteristics

Type of collection	venous
Type of preparation	Clotting activator/gel
Colour code	EU/ISO
Label/ Print	with plastic label
Colour of print/label	brown
Graduation	yes
Cap	flat cap
Closure type	push cap
Application	venous blood collection
Base shape	skirted conical base
Volume of work	1,1 ml

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Micro sample tube Serum Gel CAT, 1.1 ml, push cap, EU/ISO

Page 2

Recommended centrifugation 10000 U/min
Centrifugation time 5 Min

Size

Sample volume 1.1 ml
Diameter 10.8 mm
Length of product 44 mm
Length including cap 44 mm
Length excluding cap 40 mm

Material & colours

Colour of product transparent
Closure material Low Density Polyethylene (LD-PE)
Colour of cap brown
Tube material Polypropylene (PP)

Purity & certification

Product category In-vitro diagnostic (IVD)
CE certified CE - manufacturer's self-declaration
Batched yes

Packaging

Minimum order qty. 1000
Type of smallest subpackaging bag
Piece(s) / inner box 100
Piece (s) / outer case 1000
Piece(s) / pallet 96000
Depth of case 348 mm
Width of case 278 mm
Height of case 140 mm
Case volume 0.0135 cbm
Weight of product 0.0018 kg
Weight of case 2.11 kg
EAN of inner box 4038917129473
EAN case 4038917022248

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

