



**Product Specification**  
**S-Monovette® Serum CAT, 7.5 ml, cap white,**  
**(LxØ): 92 x 15 mm, with paper label**

**Page 1**



**Product description**

Order number	01.1601.001
Product description	S-Monovette® Serum CAT, preparation: Clotting activator, 7.5 ml, membrane screw cap, cap white, colour code EU, (LxØ) without cap: 92 x 15 mm, with paper label, label/print: white, 50 piece(s)/case, sterile

**Product characteristics**

Type of collection	venous
Type of preparation	Clotting activator
Colour code	EU
Label/ Print	with paper label
Colour of print/label	white
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2000 x g
Centrifugation temperature	20
Centrifugation time	10 Min
Liquid-tight	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.







**Product Specification**  
**S-Monovette® Serum CAT, 7.5 ml, cap white,**  
**(LxØ): 92 x 15 mm, with paper label**

**Page 3**

Depth of case	400 mm
Width of case	144 mm
Height of case	440 mm
Case volume	0.0253 cbm
Weight of product	0.0084 kg
Weight of case	5.04 kg
EAN of inner box	4038917131575
EAN case	4038917023788

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-30

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**