



Product Specification
S-Monovette® Serum CAT, 7.5 ml, cap black,
(LxØ): 92 x 15 mm, with paper label

Page 1



Product description

| | |
|---------------------|--|
| Order number | 01.1601.901 |
| Product description | S-Monovette® Serum CAT, preparation: Clotting activator, 7.5 ml, membrane screw cap, cap black, (LxØ) without cap: 92 x 15 mm, with paper label, label/print: black, 50 piece(s)/case, sterile |

Product characteristics

| | |
|----------------------------|--------------------|
| Type of collection | venous |
| Type of preparation | Clotting activator |
| Label/ Print | with paper label |
| Colour of print/label | black |
| Graduation | fill mark |
| Cap | membrane screw cap |
| Closure type | screw cap |
| Recommended centrifugation | 2000 x g |
| Centrifugation temperature | 20 |
| Centrifugation time | 10 Min |
| Liquid-tight | yes |
| Connection | membrane screw cap |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
S-Monovette® Serum CAT, 7.5 ml, cap black,
(LxØ): 92 x 15 mm, with paper label

Page 2

Size

| | |
|----------------------------------|--------|
| Sample volume | 7.5 ml |
| Diameter | 15 mm |
| Length including cap | 108 mm |
| Length excluding cap | 92 mm |
| Length including cap and plunger | 129 mm |

Material & colours

| | |
|-------------------------|-----------------------------------|
| Colour of product | transparent |
| Closure material | High Density Polyethylene (HD-PE) |
| Colour of cap | black |
| Membrane material | rubber |
| Tube material | Polypropylene (PP) |
| Piston material | High Density Polyethylene (HD-PE) |
| Piston plunger material | Polystyrene (PS) |

Purity & certification

| | |
|---------------------------|--------------------------------|
| Satisfies the requirement | IATA, ADR |
| Product category | In vitro diagnostic device, CE |
| CE certified | CE |
| Purity standard | sterile |
| Sterilisation | Electron irradiation |
| Batched | yes |

Packaging

| | |
|-------------------------------|--------|
| Minimum order qty. | 500 |
| Type of smallest subpackaging | case |
| Piece(s) / inner box | 50 |
| Piece (s) / outer case | 500 |
| Piece(s) / pallet | 32000 |
| Depth of box | 189 mm |
| Width of box | 134 mm |
| Height of box | 82 mm |
| Depth of case | 400 mm |
| Width of case | 144 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
S-Monovette® Serum CAT, 7.5 ml, cap black,
(LxØ): 92 x 15 mm, with paper label

Page 3

| | |
|-------------------|---------------|
| Height of case | 440 mm |
| Case volume | 0.0253 cbm |
| Weight of product | 0.0084 kg |
| Weight of case | 5.05 kg |
| EAN of inner box | 4038917131780 |
| EAN case | 4038917023962 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-29

This document was prepared by EDP support and is valid without signature.



SARSTEDT