



Product Specification
S-Monovette® Serum Gel CAT, 7.5 ml, cap brown,
(LxØ): 92 x 15 mm, with paper label

Page 1



Product description

Order number	01.1602.001
Product description	S-Monovette® Serum Gel CAT, preparation: Clotting activator/gel, 7.5 ml, membrane screw cap, cap brown, colour code EU/ISO, (LxØ) without cap: 92 x 15 mm, with paper label, label/print: brown, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	Clotting activator/gel
Colour code	EU/ISO
Label/ Print	with paper label
Colour of print/label	brown
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2500 x g
Centrifugation temperature	20
Centrifugation time	10 Min
Liquid-tight	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
S-Monovette® Serum Gel CAT, 7.5 ml, cap brown,
(LxØ): 92 x 15 mm, with paper label

Page 2

Connection membrane screw cap

Size

Sample volume	7.5 ml
Outer diameter of cap	18 mm
Diameter	15 mm
Length including cap	108 mm
Length excluding cap	92 mm
Length including cap and plunger	129 mm

Material & colours

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	brown
Membrane material	rubber
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order qty.	500
Type of smallest subpackaging	case
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	24000
Depth of box	189 mm
Width of box	134 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
S-Monovette® Serum Gel CAT, 7.5 ml, cap brown,
(LxØ): 92 x 15 mm, with paper label

Page 3

Height of box	82 mm
Depth of case	435 mm
Width of case	400 mm
Height of case	146 mm
Case volume	0.0254 cbm
Weight of product	0.0088 kg
Weight of case	5.24 kg
EAN of inner box	4038917131797
EAN case	4038917023979

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-12-05

This document was prepared by EDP support and is valid without signature.



SARSTEDT