



Product description

Order number	01.1604
Product description	S-Monovette® Lithium heparin LH, preparation: Lithium heparin, 7.5 ml, membrane screw cap, cap orange, colour code EU, (LxØ) without cap: 92 x 15 mm, with plastic label, label/print: transparent/orange, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	Lithium heparin
Preparation concentration	16 IU/ml blood
Colour code	EU
Label/ Print	with plastic label
Colour of print/label	transparent/orange
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2000 x g
Centrifugation temperature	20
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
S-Monovette® Lithium heparin LH, 7.5 ml, cap orange, (LxØ): 92 x 15 mm, with plastic label

Page 2

Liquid-tight yes
Connection membrane screw cap

Size

Sample volume 7.5 ml
Diameter 15 mm
Length including cap 108 mm
Length excluding cap 92 mm
Length including cap and plunger 129 mm

Material & colours

Colour of product transparent
Closure material High Density Polyethylene (HD-PE)
Colour of cap orange
Membrane material rubber
Tube material Polypropylene (PP)
Piston material High Density Polyethylene (HD-PE)
Piston plunger material Polystyrene (PS)

Purity & certification

Satisfies the requirement IATA, ADR
Product category In vitrodiagnostic device, CE
CE certified CE
Purity standard sterile
Sterilisation Electron irradiation
Batched yes

Packaging

Minimum order qty. 500
Type of smallest case
subpackaging
Piece(s) / inner box 50
Piece (s) / outer case 500
Piece(s) / pallet 32000
Depth of box 189 mm
Width of box 134 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
S-Monovette® Lithium heparin LH, 7.5 ml, cap orange, (LxØ): 92 x 15 mm, with plastic label

Height of box	82 mm
Depth of case	440 mm
Width of case	400 mm
Height of case	129 mm
Case volume	0.0227 cbm
Weight of product	0.0084 kg
Weight of case	5.03 kg
EAN of inner box	4038917087223
EAN case	4038917005692

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

