



## Product Specification

**S-Monovette® Ammonium heparin AH, 9 ml, cap blue, (LxØ): 92 x 16 mm, with plastic label**

**Page 1**



### Product description

Order number	02.1064
Product description	S-Monovette® Ammonium heparin AH, preparation: Ammonium heparin, 9 ml, membrane screw cap, cap blue, colour code EU, (LxØ) without cap: 92 x 16 mm, with plastic label, label/print: transparent/blue, 50 piece(s)/case, sterile

### Product characteristics

Type of collection	venous
Type of preparation	Ammonium heparin
Preparation concentration	16 IU/ml blood
Colour code	EU
Label/ Print	with plastic label
Colour of print/label	transparent/blue
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2000 x g
Centrifugation temperature	20
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-05

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**S-Monovette® Ammonium heparin AH, 9 ml, cap blue, (LxØ): 92 x 16 mm, with plastic label**

**Page 2**

Liquid-tight Connection	yes membrane screw cap
-------------------------	---------------------------

### Size

Sample volume	9 ml
Diameter	16 mm
Length including cap	108 mm
Length excluding cap	92 mm
Length including cap and plunger	129 mm

### Material & colours

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	blue
Membrane material	rubber
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

### Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

### Packaging

Minimum order qty.	500
Type of smallest subpackaging	case
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	24000
Depth of box	204 mm
Width of box	134 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-05

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**

**S-Monovette® Ammonium heparin AH, 9 ml, cap blue, (LxØ): 92 x 16 mm, with plastic label**

Height of box	88 mm
Depth of case	648 mm
Width of case	303 mm
Height of case	150 mm
Case volume	0.0295 cbm
Weight of product	0.0088 kg
Weight of case	5.33 kg
EAN of inner box	4038917085229
EAN case	4038917003919

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

