



**Product Specification**  
**S-Monovette® Serum Gel CAT, 4 ml, cap brown,**  
**(LxØ): 75 x 13 mm, with paper label**

**Page 1**



**Product description**

Order number	04.1925.001
Product description	S-Monovette® Serum Gel CAT, preparation: Clotting activator/gel, 4 ml, membrane screw cap, cap brown, colour code EU/ISO, (LxØ) without cap: 75 x 13 mm, with paper label, label/print: brown, 50 piece(s)/case, sterile

**Product characteristics**

Type of collection	venous
Type of preparation	Clotting activator/gel
Colour code	EU/ISO
Label/ Print	with paper label
Colour of print/label	brown
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2500 x g
Centrifugation temperature	20
Centrifugation time	10 Min
Liquid-tight	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.







**Product Specification**  
**S-Monovette® Serum Gel CAT, 4 ml, cap brown,**  
**(LxØ): 75 x 13 mm, with paper label**

Height of box	73 mm
Depth of case	388 mm
Width of case	358 mm
Height of case	135 mm
Case volume	0.0188 cbm
Weight of product	0.0055 kg
Weight of case	3.49 kg
EAN of inner box	4038917359443
EAN case	4038917359436

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

