



**Product Specification**  
**S-Monovette® EDTA K3E, 4.9 ml, cap red, (LxØ):**  
**90 x 13 mm, with paper label**

**Page 1**



**Product description**

Order number	04.1931.001
Product description	S-Monovette® EDTA K3E, preparation: K3 EDTA, 4.9 ml, membrane screw cap, cap red, colour code EU, (LxØ) without cap: 90 x 13 mm, with paper label, label/print: red, 50 piece(s)/case, sterile

**Product characteristics**

Type of collection	venous
Type of preparation	K3 EDTA
Preparation concentration	1.6 mg/ml blood
Colour code	EU
Label/ Print	with paper label
Colour of print/label	red
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Liquid-tight	yes
Connection	membrane screw cap

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**S-Monovette® EDTA K3E, 4.9 ml, cap red, (LxØ):**  
**90 x 13 mm, with paper label**

**Page 2**

**Size**

Sample volume	4.9 ml
Diameter	13 mm
Length including cap	106 mm
Length excluding cap	90 mm
Length including cap and plunger	127 mm

**Material & colours**

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	red
Membrane material	rubber
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

**Purity & certification**

Satisfies the requirement	IATA, ADR
Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

**Packaging**

Minimum order qty.	500
Type of smallest subpackaging	case
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	36000
Depth of box	171 mm
Width of box	130 mm
Height of box	73 mm
Depth of case	388 mm
Width of case	358 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**S-Monovette® EDTA K3E, 4.9 ml, cap red, (LxØ):**  
**90 x 13 mm, with paper label**

Height of case	150 mm
Case volume	0.0208 cbm
Weight of product	0.0049 kg
Weight of case	3.21 kg
EAN of inner box	4038917134446
EAN case	4038917025690

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

