



Product description

Order number	04.2911.200
Product description	S-Monovette® Serum CAT, CAT, preparation: Clotting activator, 1.6 ml, membrane screw cap, cap red, colour code ISO, (LxØ) without cap: 75 x 13 mm, with plastic label, label/print: transparent/red, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	Clotting activator
Colour code	ISO
Label/ Print	with plastic label
Colour of print/label	transparent/red
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2000 x g
Centrifugation temperature	20
Centrifugation time	10 Min
Liquid-tight	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
S-Monovette® Serum CAT, 1.6 ml, cap red, (LxØ):
75 x 13 mm, with plastic label

Page 2

Connection membrane screw cap

Size

Sample volume	1.6 ml
Outer diameter of cap	16 mm
Diameter	13 mm
Length including cap	91 mm
Length excluding cap	75 mm
Length including cap and plunger	112 mm

Material & colours

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	red
Membrane material	rubber
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order qty.	500
Type of smallest subpackaging	case
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	39000
Depth of box	171 mm
Width of box	116 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
S-Monovette® Serum CAT, 1.6 ml, cap red, (LxØ):
75 x 13 mm, with plastic label

Height of box	73 mm
Depth of case	388 mm
Width of case	358 mm
Height of case	135 mm
Case volume	0.0188 cbm
Weight of product	0.0049 kg
Weight of case	3.18 kg
EAN of inner box	4038917983327
EAN case	4038917983334

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

