



Product description

| | |
|---------------------|---|
| Order number | 05.1557.100 |
| Product description | S-Monovette® Serum CAT, preparation: Clotting activator, 2.7 ml, membrane screw cap, cap red, colour code ISO, (LxØ) without cap: 66 x 11 mm, with paper label, label/print: red, 50 piece(s)/case, sterile |

Product characteristics

| | |
|----------------------------|--------------------|
| Type of collection | venous |
| Type of preparation | Clotting activator |
| Colour code | ISO |
| Label/ Print | with paper label |
| Colour of print/label | red |
| Graduation | fill mark |
| Cap | membrane screw cap |
| Closure type | screw cap |
| Recommended centrifugation | 2000 x g |
| Centrifugation temperature | 20 |
| Centrifugation time | 10 Min |
| Liquid-tight | yes |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
S-Monovette® Serum CAT, 2.7 ml, cap red, (LxØ):
66 x 11 mm, with paper label

Page 2

Connection membrane screw cap

Size

| | |
|----------------------------------|--------|
| Sample volume | 2.7 ml |
| Outer diameter of cap | 13 mm |
| Diameter | 11 mm |
| Length including cap | 82 mm |
| Length excluding cap | 66 mm |
| Length including cap and plunger | 103 mm |

Material & colours

| | |
|-------------------------|-----------------------------------|
| Colour of product | transparent |
| Closure material | High Density Polyethylene (HD-PE) |
| Colour of cap | red |
| Membrane material | rubber |
| Tube material | Polypropylene (PP) |
| Piston material | High Density Polyethylene (HD-PE) |
| Piston plunger material | Polystyrene (PS) |

Purity & certification

| | |
|---------------------------|--------------------------------------|
| Satisfies the requirement | IATA, ADR |
| Product category | In-vitro diagnostic (IVD) |
| CE certified | CE - manufacturer's self-declaration |
| Purity standard | sterile |
| Sterilisation | Electron irradiation |
| Batched | yes |

Packaging

| | |
|-------------------------------|--------|
| Minimum order qty. | 500 |
| Type of smallest subpackaging | case |
| Piece(s) / inner box | 50 |
| Piece (s) / outer case | 500 |
| Piece(s) / pallet | 56000 |
| Depth of box | 145 mm |
| Width of box | 108 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
S-Monovette® Serum CAT, 2.7 ml, cap red, (LxØ):
66 x 11 mm, with paper label

| | |
|-------------------|---------------|
| Height of box | 63 mm |
| Depth of case | 333 mm |
| Width of case | 303 mm |
| Height of case | 125 mm |
| Case volume | 0.0126 cbm |
| Weight of product | 0.0042 kg |
| Weight of case | 2.63 kg |
| EAN of inner box | 4038917130691 |
| EAN case | 4038917032315 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

