




Product Specification
S-Monovette® Fluoride/EDTA FE, 1.2 ml, cap yellow, (LxØ): 66 x 8 mm, with plastic label

Page 1



Product description

Order number	06.1665.001
Product description	S-Monovette® Fluoride/EDTA FE, preparation: Fluoride + EDTA, 1.2 ml, membrane screw cap, cap yellow, colour code EU, (LxØ) without cap: 66 x 8 mm, with plastic label, label/print: yellow/black, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	Fluoride + EDTA
Preparation concentration	1.0 mg/ml blood + 1.2 mg/ml blood
Colour code	EU
Label/ Print	with plastic label
Colour of print/label	yellow/black
Graduation	yes
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2000 x g
Centrifugation temperature	20
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

S-Monovette® Fluoride/EDTA FE, 1.2 ml, cap yellow, (LxØ): 66 x 8 mm, with plastic label

Page 2

Liquid-tight Connection	yes membrane screw cap
-------------------------	---------------------------

Size

Sample volume	1.2 ml
Outer diameter of cap	13 mm
Diameter	8 mm
Length including cap	82 mm
Length excluding cap	66 mm
Length including cap and plunger	103 mm

Material & colours

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	yellow
Membrane material	rubber
Tube material	Polypropylene (PP)
Piston material	rubber
Piston plunger material	Polystyrene (PS)

Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order qty.	500
Type of smallest subpackaging	case
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	56000
Depth of box	145 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-12-10

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
**S-Monovette® Fluoride/EDTA FE, 1.2 ml, cap
yellow, (LxØ): 66 x 8 mm, with plastic label**

Page 3

Width of box	108 mm
Height of box	63 mm
Depth of case	333 mm
Width of case	303 mm
Height of case	125 mm
Case volume	0.0126 cbm
Weight of product	0.0033 kg
Weight of case	2.19 kg
EAN of inner box	4038917132848
EAN case	4038917024488

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-12-10

This document was prepared by EDP support and is valid without signature.



SARSTEDT