



## Product Specification

**S-Monovette® Citrate 9NC 0.106 mol/l 3.2%, 1.4 ml, cap blue, (LxØ): 66 x 8 mm, with plastic label**

**Page 1**



### Product description

Order number	06.1668.100
Product description	S-Monovette® Citrate 9NC 0.106 mol/l 3.2%, preparation: Citrate 9NC, 1.4 ml, membrane screw cap, cap blue, colour code ISO, (LxØ) without cap: 66 x 8 mm, with plastic label, label/print: white/blue, 50 piece(s)/case, sterile

### Product characteristics

Type of collection	venous
Type of preparation	Citrate 9NC
Preparation concentration	0.106 mol/l
Colour code	ISO
Label/ Print	with plastic label
Colour of print/label	white/blue
Graduation	yes
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	1800 x g
Centrifugation temperature	20
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-10

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**S-Monovette® Citrate 9NC 0.106 mol/l 3.2%, 1.4 ml, cap blue, (LxØ): 66 x 8 mm, with plastic label**

**Page 2**

Liquid-tight Connection      yes  
membrane screw cap

### Size

Sample volume      1.4 ml  
Diameter      8 mm  
Length including cap      82 mm  
Length excluding cap      66 mm  
Length including cap and plunger      103 mm

### Material & colours

Colour of product      transparent  
Closure material      High Density Polyethylene (HD-PE)  
Colour of cap      blue  
Membrane material      rubber  
Tube material      Polypropylene (PP)  
Piston material      rubber  
Piston plunger material      Polystyrene (PS)

### Purity & certification

Satisfies the requirement      IATA, ADR  
Product category      In vitrodiagnostic device, CE  
CE certified      CE  
Purity standard      sterile  
Sterilisation      Electron irradiation  
Batched      yes

### Packaging

Minimum order qty.      500  
Type of smallest subpackaging      case  
Piece(s) / inner box      50  
Piece (s) / outer case      500  
Piece(s) / pallet      56000  
Depth of box      145 mm  
Width of box      108 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-10

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**

**S-Monovette® Citrate 9NC 0.106 mol/l 3.2%, 1.4 ml, cap blue, (LxØ): 66 x 8 mm, with plastic label**

Height of box	63 mm
Depth of case	333 mm
Width of case	303 mm
Height of case	125 mm
Case volume	0.0126 cbm
Weight of product	0.0033 kg
Weight of case	2.19 kg
EAN of inner box	4038917132916
EAN case	4038917024556

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

