



## Product description

Order number	04.1929
Product description	S-Monovette® Lithium heparin LH, preparation: Lithium heparin LH, 2.7 ml, membrane screw cap, cap orange, colour code EU, (LxØ) without cap: 75 x 13 mm, with plastic label, label/print: transparent/orange, 50 piece(s)/case, sterile

## Product characteristics

Type of collection	venous
Type of preparation	Lithium heparin LH
Preparation concentration	16 IU/ml blood
Colour code	EU
Label/ Print	with plastic label
Colour of print/label	transparent/orange
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2000 x g
Centrifugation temperature	20

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**S-Monovette® Lithium heparin LH, 2.7 ml, cap orange, (LxØ): 75 x 13 mm, with plastic label**

**Page 2**

Centrifugation time	10 Min
Liquid-tight	yes
Connection	membrane screw cap

**Size**

Sample volume	2.7 ml
Diameter	13 mm
Length including cap	91 mm
Length excluding cap	75 mm
Length including cap and plunger	112 mm

**Material & colours**

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	orange
Membrane material	rubber
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

**Purity & certification**

Satisfies the requirement	IATA, ADR
Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

**Packaging**

Minimum order qty.	500
Type of smallest subpackaging	case
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	39000
Depth of box	171 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**S-Monovette® Lithium heparin LH, 2.7 ml, cap orange, (LxØ): 75 x 13 mm, with plastic label**

Width of box	116 mm
Height of box	73 mm
Depth of case	388 mm
Width of case	358 mm
Height of case	135 mm
Case volume	0.0188 cbm
Weight of product	0.0045 kg
Weight of case	3 kg
EAN of inner box	4038917360678
EAN case	4038917360661

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

