



Product Specification

S-Monovette® Lithium heparin gel+ LH, 4 ml, cap orange, (LxØ): 75 x 13 mm, with plastic label

Page 1



Product description

| | |
|---------------------|---|
| Order number | 04.1953 |
| Product description | S-Monovette® Lithium heparin gel+ LH, preparation: Lithium heparin gel+, 4 ml, membrane screw cap, cap orange, colour code EU, (LxØ) without cap: 75 x 13 mm, with plastic label, label/print: transparent/brown, 50 piece(s)/case, sterile |

Product characteristics

| | |
|----------------------------|----------------------|
| Type of collection | venous |
| Type of preparation | Lithium heparin gel+ |
| Preparation concentration | 25 IU/ml blood |
| Colour code | EU |
| Label/ Print | with plastic label |
| Colour of print/label | transparent/brown |
| Graduation | fill mark |
| Cap | membrane screw cap |
| Closure type | screw cap |
| Recommended centrifugation | 2500 x g |
| Centrifugation temperature | 20 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-09-20

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

S-Monovette® Lithium heparin gel+ LH, 4 ml, cap orange, (LxØ): 75 x 13 mm, with plastic label

Page 2

| | |
|---------------------|--------------------|
| Centrifugation time | 7 Min |
| Liquid-tight | yes |
| Connection | membrane screw cap |

Size

| | |
|----------------------------------|--------|
| Sample volume | 4 ml |
| Outer diameter of cap | 16 mm |
| Diameter | 13 mm |
| Length including cap | 91 mm |
| Length excluding cap | 75 mm |
| Length including cap and plunger | 112 mm |

Material & colours

| | |
|-------------------------|-----------------------------------|
| Colour of product | transparent |
| Closure material | High Density Polyethylene (HD-PE) |
| Colour of cap | orange |
| Membrane material | rubber |
| Tube material | Polypropylene (PP) |
| Piston material | High Density Polyethylene (HD-PE) |
| Piston plunger material | Polystyrene (PS) |

Purity & certification

| | |
|---------------------------|--------------------------------|
| Satisfies the requirement | IATA, ADR |
| Product category | In vitro diagnostic device, CE |
| CE certified | CE |
| Purity standard | sterile |
| Sterilisation | Electron irradiation |
| Batched | yes |

Packaging

| | |
|-------------------------------|-------|
| Minimum order qty. | 500 |
| Type of smallest subpackaging | case |
| Piece(s) / inner box | 50 |
| Piece (s) / outer case | 500 |
| Piece(s) / pallet | 39000 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-09-20

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

S-Monovette® Lithium heparin gel+ LH, 4 ml, cap orange, (LxØ): 75 x 13 mm, with plastic label

Page 3

| | |
|-------------------|---------------|
| Depth of box | 171 mm |
| Width of box | 116 mm |
| Height of box | 73 mm |
| Depth of case | 388 mm |
| Width of case | 358 mm |
| Height of case | 135 mm |
| Case volume | 0.0188 cbm |
| Weight of product | 0.0054 kg |
| Weight of case | 3.45 kg |
| EAN of inner box | 4038917424929 |
| EAN case | 4038917424912 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-09-20

This document was prepared by EDP support and is valid without signature.



SARSTEDT