




### **Product description**

|                     |                                                                                                                 |
|---------------------|-----------------------------------------------------------------------------------------------------------------|
| Order number        | 83.1840.001                                                                                                     |
| Product description | Coverslips, for Adherent cells, Ø: 22 mm, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 25 piece(s)/bag |

### **Product characteristics**

|                           |                |
|---------------------------|----------------|
| Suitable for              | 6-well plates  |
| Cell type                 | Adherent cells |
| Type of surface treatment | TC-treated     |
| Coding colour             | red            |

### **Size**

|          |       |
|----------|-------|
| Diameter | 22 mm |
|----------|-------|

### **Material & colours**

|                   |                                                     |
|-------------------|-----------------------------------------------------|
| Product material  | Polyethylene terephthalate, Glycol-modified (PET-G) |
| Colour of product | transparent                                         |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

|                  |                                                      |
|------------------|------------------------------------------------------|
| Product category | no medical device   IVD                              |
| Purity standard  | sterile, non-pyrogenic/endotoxin-free, non-cytotoxic |
| Sterilisation    | Gamma irradiation                                    |
| Pyrogen-free     | yes                                                  |
| Non-cytotoxic    | yes                                                  |
| Batched          | yes                                                  |

Packaging

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 200           |
| Type of smallest subpackaging | bag           |
| Piece(s) / inner box          | 200           |
| Piece (s) / outer case        | 2000          |
| Piece(s) / pallet             | 72000         |
| Depth of case                 | 310 mm        |
| Width of case                 | 144 mm        |
| Height of case                | 335 mm        |
| Case volume                   | 0.015 cbm     |
| Weight of product             | 0.0003 kg     |
| Weight of case                | 1.8 kg        |
| EAN of inner box              | 4038917134040 |
| EAN case                      | 4038917025386 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

