



Product description

Order number	83.1840.002
Product description	Coverslips, for Adherent cells, Ø: 13 mm, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 25 piece(s)/bag

Product characteristics

Suitable for	6-, 12- and 24-well plates
Cell type	Adherent cells
Type of surface treatment	TC-treated
Coding colour	red

Size

Diameter	13 mm
----------	-------

Material & colours

Product material	Polyethylene terephthalate, Glycol-modified (PET-G)
Colour of product	transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Product category	no medical device IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

Packaging

Minimum order qty.	200
Type of smallest subpackaging	bag
Piece(s) / inner box	200
Piece (s) / outer case	2000
Piece(s) / pallet	288000
Depth of case	310 mm
Width of case	144 mm
Height of case	335 mm
Case volume	0.015 cbm
Weight of product	0.0003 kg
Weight of case	1.8 kg
EAN of inner box	4038917134057
EAN case	4038917025393

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

