




Product Specification
Serological pipette, with tip, plugged, 5 ml, sterile,
25 piece(s)/bag

Page 1



Product description

Order number	86.1253.025
Product description	Serological pipette, with tip, plugged, 5 ml, coding colour: blue, material: PS, sterile, 25 piece(s)/bag

Product characteristics

Colour of print/label	blue
Graduation	yes
Design	with tip, plugged
Volume of work	5 ml
Coding colour	blue

Size

Diameter	8 mm
Length of product	340 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent

Purity & certification

Product category	no medical device IVD
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes

Packaging

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	25
Piece (s) / outer case	500
Piece(s) / pallet	40500
Depth of case	413 mm
Width of case	250 mm
Height of case	195 mm
Case volume	0.0201 cbm
Weight of product	0.0067 kg
Weight of case	4.03 kg
EAN of inner box	4038917127257
EAN case	4038917021111

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

