




**Product Specification**  
**Serological pipette, with tip, plugged, 10 ml,  
sterile, non-pyrogenic/endotoxin-free,  
non-cytotoxic, 1 piece(s)/blister**

**Page 1**



**Product description**

|                     |                                                                                                                                                              |
|---------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Order number        | 86.1254.001                                                                                                                                                  |
| Product description | Serological pipette, with tip, plugged, 10 ml, coding colour: orange, material: PS, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister |

**Product characteristics**

|                       |                   |
|-----------------------|-------------------|
| Colour of print/label | orange            |
| Graduation            | yes               |
| Design                | with tip, plugged |
| Volume of work        | 10 ml             |
| Coding colour         | orange            |

**Size**

|                   |        |
|-------------------|--------|
| Diameter          | 10 mm  |
| Length of product | 335 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Serological pipette, with tip, plugged, 10 ml,  
sterile, non-pyrogenic/endotoxin-free,  
non-cytotoxic, 1 piece(s)/blister**

**Page 2**

**Material & colours**

|                   |                  |
|-------------------|------------------|
| Product material  | Polystyrene (PS) |
| Colour of product | transparent      |

**Purity & certification**

|                  |                                                      |
|------------------|------------------------------------------------------|
| Product category | no medical device   IVD                              |
| Purity standard  | sterile, non-pyrogenic/endotoxin-free, non-cytotoxic |
| Sterilisation    | Gamma irradiation                                    |
| Pyrogen-free     | yes                                                  |
| Batched          | yes                                                  |

**Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 500           |
| Type of smallest subpackaging | blister       |
| Piece(s) / inner box          | 50            |
| Piece (s) / outer case        | 500           |
| Piece(s) / pallet             | 27000         |
| Depth of case                 | 420 mm        |
| Width of case                 | 250 mm        |
| Height of case                | 290 mm        |
| Case volume                   | 0.0305 cbm    |
| Weight of product             | 0.0104 kg     |
| Weight of case                | 5.74 kg       |
| EAN of inner box              | 4038917127295 |
| EAN case                      | 4038917021128 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-29

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**