



Product description

Order number	94.6077.435
Product description	flexiPERM® , con B, re-usable silicon insert, for dividing the growth surface, 1 well, non-pyrogenic/endotoxin-free, non-cytotoxic, 5 piece(s)/bag

Product characteristics

Application	for dividing the growth surface
Number of wells	1
Growth area	3.1 cm ²
Volume of work	2 ml - 3 ml

Size

Diameter	20 mm
Height of product	8 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	silicone
Colour of product	transparent

Purity & certification

Product category	no medical device IVD
Purity standard	non-pyrogenic/endotoxin-free, non-cytotoxic
Non-cytotoxic	yes
Batched	yes

Packaging

Minimum order qty.	5
Type of smallest subpackaging	bag
Piece(s) / inner box	5
Piece (s) / outer case	100
Piece(s) / pallet	7200
Depth of case	388 mm
Width of case	253 mm
Height of case	225 mm
Case volume	0.0221 cbm
Weight of product	0.0039 kg
Weight of case	1.42 kg
EAN of inner box	4038917255295
EAN case	4038917255288

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

