



### Product description

Order number	83.3925.400
Product description	Cell culture plate, 96 well, material: PS, surface: BIOFLOAT?, for spheroid culture, round base, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/bag

### Product characteristics

Closure type	slip-on lid
Cell type	spheroid culture
Base shape	round base
Number of wells	96
Volume of work	0,1 ml - 0,2 ml
Vessel type	Plate

### Size

Max. volume	0.31 ml
Width of product	85.2 mm
Height of product	16.55 mm
Length of product	127.8 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent
Closure material	Polystyrene (PS)

**Purity & certification**

Product category	no medical device   IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Electron irradiation
Non-cytotoxic	yes
Batched	yes

**Packaging**

Minimum order qty.	4
Type of smallest subpackaging	bag
Piece(s) / inner box	4
Piece (s) / outer case	24
Piece(s) / pallet	1440
Depth of case	322 mm
Width of case	257 mm
Height of case	245 mm
Case volume	0.0203 cbm
Weight of product	0.049 kg
Weight of case	2.36 kg
EAN of inner box	4038917923965
EAN case	4038917923972

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

